

ISO GUIDE 34

ISO/IEC 17025

ISO 9001:2008

Certificate of Analysis

Equilin

Catalog Number:	Е-078 нас и
Solution Lot:	FN080411-02
Retest Date:	November 2014
Solvent:	Methanol
Volume per Ampule:	Not less than 1 mL
Storage:	Store in Freezer. See Stability Section.
Intended Use:	For R&D/ analytical purposes only. Not suitable for human or animal consumption.
Safety:	Flammable, Poison

Retest Date - stability studies ongoing. Certificate of Analysis will be updated upon completion of retest.

• Ampules are overfilled to ensure a minimum 1 mL volume fill. We advise laboratories to use measured volumes of this standard solution before diluting to the desired concentration.

	Component	Chromatographic Purity	Certified Concentration
	Equilin	99.9%	1.000 ± 0.006 mg/mL
٠	Uncertainty of the concentration is ex	pressed as an expanded uncertainty in accordance with	ISO 17025 and Guide 34 at the approximate
	95% confidence interval using a coverage factor of k = 2 and has been calculated by statistical analysis of our production system and		
	incorporates uncertainty of the purity factor, material density, and balance and weighing technique.		
٠	Concentration is corrected for chroma	tographic purity, residual water, residual solvents and	residual inorganics.

Solution Standard Verification and Homogeneity

Standard		Verified Concentration (mg/mL)		%RSD - Homogeneity	
Solution	Lot Number	Actual Results	Acceptance Criteria	Actual Results	Acceptance Criteria
New Lot	FN080411-02	1.003	± 3%	0.8	≤ 3%
• Concentration is verified through multiple analyses and is calculated as the average of multiple analyses compared to an independently					

Concentration is verified through multiple analyses and is calculated as the average of multiple analyses compared to an independently
prepared calibration solution.

 Homogeneity of the New Lot is ensured through rigorous production process controls statistically analyzed to evaluate risk and verified by analysis. The % RSD of samples pulled from across the lot demonstrate homogeneity of the New Lot.

Traceability

- Gravimetrically prepared using qualified balances calibrated semi-annually by Mettler Toledo using NIST traceable weights. Calibration
 verification performed weekly and prior to each use utilizing NIST traceable weights. Each balance has been assigned a minimum weighing
 by Mettler Toledo taking into consideration the balance and installed environmental conditions to ensure weighing complies with USP
 tolerances of no more than 0.1% relative error.
- Concentration is verified against an independently prepared calibration solution gravimetrically prepared using balances calibrated to NIST.
- In addition, each neat material utilized has been identified and thoroughly characterized through the use of multiple analytical techniques. Spectral data is provided on subsequent pages of the COA.

Cerilliant certifies that this standard meets the specifications stated in this certificate and warrants this product to meet the stated acceptance criteria through the expiration/retest date when stored unopened as recommended. Product should be used shortly after opening to avoid concentration changes due to evaporation. Warranty does not apply to ampoules stored after opening.



Darron Ellsworth, Quality Assurance Manager

October 15, 2013

Date

E-078 FN080411-02 Revision 2 Page 2 of 6

99.84%

Standard Soluti	ion Assay Parameters	Calibration (Calibration Curve		
Analysis Method: Wavelength: Slit Width: Response:	UV/Vis 281 nm 1.0 nm 0.5 s	Calibration Curve: Number of Points: Linearity (r):	Linear Regression 4 1.000		
Neat Material L	Data				
Compound Name: Compound Lot:	Equilin PN071111-01	Chemical Formula: CAS Number: Molecular Weight:	C ₁₈ H ₂₀ O ₂ 474-86-2 268.35		
	Neat Material C	Characterization Summary			
Analytical Test		Method	Results		
Primary Chromatographi	c Purity by HPLC/PDA Analysis	SP10-0102	99.9%		
Secondary Chromatogra	phic Purity by GC/FID Analysis	SP10-0101	99.6%		
dentity by GC/MS Anal	ysis	SP10-0105	Consistent with Structure		
dentity by ¹ H-NMR Ana	alysis	USP <761>, SP10-0116	Consistent with Structure		
Residual Solvent Analysi	is by GC/FID Headspace	AM1087 ¹	None Detected		
Residual Water Analysis	by Karl Fischer Coulometry	USP <921>, SP10-0103	0.07%		
norganic Content by Mi	croash Analysis	SP10-0135	< 0.2%		

• Primary purity is calculated as the average of two independently performed analyses utilizing two different methods. Acceptance criteria requires the purity values to be within 0.5% of each other.

• The primary chromatographic purity value is used to calculate the Purity Factor.

• A secondary chromatographic purity method is utilized as a control.

Purity Factor = [(100 - wt% residual solvent - wt% residual water - wt% residual inorganics) x Chromatographic Purity/100].

• Purity factor does not include adjustment for chiral and/or isotopic purity.

¹ Validated analytical method

Purity Factor

Spectral and Physical Data

HPLC/PDA



GC/FID



Column:	DB-5m	s, 30 m x 0.53 mi	m ID, 1.5 μm filı	n thickness
Temp Program:	40°C to 200°C at 40°C/min			
	200°C t	o 280°C at 5°C/n	nin (hold 18 mir	ı)
Injector Temp:	Cool-or	n-Column		
Detector Temp:	325°C	325°C		
Data File Name:	S:\GC\0	GC6\2011\0711\F	70719130.D	
Operator:	PT			
Instrument:	GC#6			
Sample Name:	PN0711	11-01		
Method File:	B01E.N	1		
Acquired:	July 20,	2011 1:49 AM		
Peak # R	et Time	Area	Height	Area %

reak #	Ret Thile	Area	neigin	Area 70
1	6.72	3.31	1.09	0.01
2	17.13	8.53	0.67	0.04
3	17.66	23872.80	4003.25	99.65
4	18.69	26.22	3.66	0.11
5	18.81	38.22	5.45	0.16
6	22.88	7.75	1.10	0.03

E-078 FN080411-02 Revision 2 Page 4 of 6

Spectral and Physical Data (cont.)





Spectral and Physical Data (cont.)

Identity by GC/MS Analysis



Stability

Short Term Stability : A summary of accelerated stability findings for the related product E-075 Estrone is listed below.

Storage Condition	Mean Kinetic Temperature (MKT)	Time Period
Freezer	-15°C	
Refrigerator	4°C	No decreases in purity was noted after one weak
Room Temperature	21°C	No decrease in purity was noted after one week.
40°C	40°C	

Transport/Shipping : Stability data supports transport of this product at MKT not exceeding 40 °C over a period of one week.

COA Revision History

Revision No.	Date	Reason for Revision
00	August 11, 2011	Initial version
01	August 20, 2012	Revised Retest Date from October 2012 to October 2013.
02	October 15, 2013	Revised Retest Date from October 2013 to November 2014.